

This document is intended to help customers to select medical flow sensors out of the Honeywell Healthcare Solutions GmbH portfolio that are compatible to the original sensors.

While Honeywell may provide information or engineering support for its products through Honeywell personnel, literature and website, it is the customer's sole responsibility to determine the suitability of the Honeywell product(s) for the customer's requirements. Specifications may change without notice. The information we supply is believed to be accurate and reliable as of this writing.

Check the Honeywell website for the latest version of this document and for technical details of the sensor specifications!

[Link to Cross Reference documents on website](https://hwll.co/envitec-cross-reference-list)

<https://hwll.co/envitec-cross-reference-list>


Information about Flow sensors

<https://automation.honeywell.com/de/de/products/sensing-solutions/sensors/flow-sensors/flow-sensors-type-a>

<https://automation.honeywell.com/de/de/products/sensing-solutions/sensors/flow-sensors/flow-sensors-type-h>

Contact: [info@envitec.com](mailto:info@envitec.com)

Original model	Original article reference	Sensor type	Honeywell part number	Honeywell sensor
Hamilton		PN 279331	E07-00-0001	SpiroQuant H 1,8 meter
Dräger	Cicero	843735 (not applying for use in Oxylog® 2000 plus)	E1000470	SpiroQuant A+
Dräger	Evita® Serie	843735 (not applying for use in Oxylog® 2000 plus)	E1000470	SpiroQuant A+
Dräger	Primus®	843735 (not applying for use in Oxylog® 2000 plus)	E1000470	SpiroQuant A+

	Revision:	Valid date:	Created:	Page – Seite:
007-35-XRL-Cross_Reference_List_Flow-0.docx	0	2025-08-22	HB	1 / 1